

HM 161.33 Crump weir



The illustration shows a similar unit

Description

■ weir according to E. S. Crump

Crump weirs are control structures. They are broad-crested weirs. The triangular shape of the weir has several advantages, e.g. only minor siltation occurs upstream of the weir. a part of sediment transport can flow over the weir. In addition, aquatic creatures can often pass this weir upstream. HM 161.33 is a weir according to E. S. Crump with defined inclinations both upstream and downstream. The Crump weir is preferably used as a sill. Sills are used to decrease the flow velocity to avoid erosion. A sill is well dimensioned for the prevailing discharge when there is no hydraulic jump.

Learning objectives/experiments

- free and submerged overfall at the Crump weir
- observation of downstream hydraulics jumps
- discharge at a sill
- together with a level gauge:
 - ► determination of the discharge
 - comparison of the theoretical and the measured discharge

Specification

- [1] Crump weir for the experimental flume HM 161
- [2] weir body contour according to E. S. Crump
- [3] weir body with sealing lips

Technical data

Weir body

- made of PVC
- inclination (upstream): 1:2
- inclination (downstream): 1:5

LxWxH: 1060x600x160mm Weight: approx. 40kg

Scope of delivery

- 1 Crump weir
- 1 set of accessories
- 1 manual

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 1/2 - 11.2023



HM 161.33 Crump weir

Required accessories

HM 161 Experimental flume 600x800mm

G.U.N.T. Gerätebau GmbH, Hanskampring 15-17, D-22885 Barsbüttel, Telefon (040) 67 08 54-0, Fax (040) 67 08 54-42, Email sales@gunt.de, Web www.gunt.de We reserve the right to modify our products without any notifications. Page 2/2 - 11.2023